

Supporting health care around the globe.



Face Shields



Face shields are used to protect the facial area and associated mucous membranes (eyes, nose, mouth) from splashes, sprays and spatter of body fluids.

PRODUCT NAME	Face Shield
PRODUCT TYPE	Fluid, droplet, large particulate protection
PRODUCT DESCRIPTION	Provide plastic barrier in front of wearers face
PERFORMANCE REQUIREMENTS	Protects facial area including eyes, nose, mouth from impact hazards and fluids
FEATURES	Washable, lightweight and adjustable, full face protection, no fit test required
SIZE	Standard - one size fits all
RECOMMENDED USE	Medical, dental, veterinary
REUSABLE	Single use recommended with exposure to contaminants
FLUID RESISTANCE	Yes
AEROSOL PROTECTION	No
LATEX FREE	Yes
CERTIFICATIONS	CE, FDA