

Supporting health care around the globe.



AAMI - Level II Isolation Gowns



Gowns are used to protect the wearer from the spread of infection or illness if in contact with potentially infectious fluids and solid material. Recommended when the expected risk of exposure to fluid is moderate.

PRODUCT NAME	AAMI - Level II Isolation Gowns
PRODUCT TYPE	Fluid and hazardous material protection
PRODUCT DESCRIPTION	Protects against hazardous particles or fluids penetration under moderate pressure conditions
PERFORMANCE REQUIREMENTS	Level II AAMI requirements met @ 4% AQL Spray impact <1.0 g, Hydrostatic pressure ≥ 20 cm
FEATURES	Layered Multi-Layer Polypropylene (SMS) fabric
SIZE	Medium, Regular (Large), X-Large, 2-X Large
RECOMMENDED USE	Hospitals - ICU, surgery, infectious control, pathology
REUSABLE	Single use
FLUID RESISTANCE	Yes - recommended for moderate exposure
AEROSOL PROTECTION	No
LATEX FREE	Yes
CERTIFICATIONS	CE, FDA, AAMI Level 2