



## AAMI - Level I Isolation Gowns



Gowns are used to protect the wearer from the spread of infection or illness if in contact with potentially infectious fluids and solid material. Recommended when the expected risk of exposure to fluid is low.

<b>PRODUCT NAME</b>	AAMI - Level I Isolation Gowns
<b>PRODUCT TYPE</b>	Fluid and hazardous material protection
<b>PRODUCT DESCRIPTION</b>	Protects against hazardous particles or fluids penetration under light pressure conditions
<b>PERFORMANCE REQUIREMENTS</b>	Level I AAMI requirements met @ 4% AQL Spray impact penetration <4.5 g
<b>FEATURES</b>	Layered Multi-Layer Polypropylene (SMS) fabric
<b>SIZE</b>	Medium, Regular (Large), X-Large, 2-X Large
<b>RECOMMENDED USE</b>	Hospitals - basic patient care
<b>REUSABLE</b>	Single use
<b>FLUID RESISTANCE</b>	Yes - recommended for moderate exposure
<b>AEROSOL PROTECTION</b>	No
<b>LATEX FREE</b>	Yes
<b>CERTIFICATIONS</b>	CE, FDA, AAMI Level 1