

Supporting health care around the globe.



N95 & KN95 Respiratory Masks



The N95 respirator is a respiratory protective device designed to achieve a very close facial fit and a 95% filtration rate of airborne particles. These surgical respirators reduce inhalation of both fluids and non-oil aerosols. The KN95 is of equivalent standard to N95 respirator masks and is a recognized alternative.

PRODUCT NAME	KN95 & N95 Respirator
PRODUCT TYPE	Aerosol filtering respirator
PRODUCT DESCRIPTION	Protects against airborne infectious agents and other airborne particles
PERFORMANCE REQUIREMENTS	95% filter efficiency level against airborne particulates $\geq 0.3 \mu\text{m}$
FEATURES	Extended wear capabilities, two-strap design for secure fit, adjustable nose clip with foam for comfort
SIZE	Standard - one size fits all
RECOMMENDED USE	Hospital - clinical use when risk of contact with airborne infections
REUSABLE	Single use recommended with exposure to contaminants
FLUID RESISTANCE	No
AEROSOL PROTECTION	Yes
LATEX FREE	Yes
CERTIFICATIONS	CE, FDA