

Supporting health care around the globe.



## NIOSH-Approved N95 Respiratory Masks



The N95 respirator is a respiratory protective device designed to achieve a very close facial fit and a 95% filtration rate of airborne particles. These surgical respirators reduce inhalation of both fluids and non-oil aerosols and are NIOSH certified with standards for air filtration efficiency, breathing resistance and leakage.

<b>PRODUCT NAME</b>	NIOSH N95 Respirator
<b>PRODUCT TYPE</b>	Aerosol filtering respirator
<b>PRODUCT DESCRIPTION</b>	Protects against airborne infectious agents and other airborne particles
<b>PERFORMANCE REQUIREMENTS</b>	NIOSH standard requirements. 95% filter efficiency level against airborne particulates >0.3 µm
<b>FEATURES</b>	Extended wear capabilities, two-strap design for secure fit, adjustable nose clip with foam for comfort
<b>SIZE</b>	Standard - one size fits all
<b>RECOMMENDED USE</b>	Hospital - clinical use when risk of contact with airborne infections
<b>REUSABLE</b>	Single use recommended with exposure to contaminants
<b>FLUID RESISTANCE</b>	No
<b>AEROSOL PROTECTION</b>	Yes
<b>LATEX FREE</b>	Yes
<b>CERTIFICATIONS</b>	CE, FDA, NIOSH