

Supporting health care around the globe.



Surgical Disposable Masks (14683 Type IIR)



Surgical masks are a barrier to splashes, droplets, and spit, and are not designed to seal tight against the face. Type IIR features fluid resistance. Specifically for medical use.

PRODUCT NAME	Surgical Masks (EN14683 - Type IIR)
PRODUCT TYPE	Fluid and droplet protection
PRODUCT DESCRIPTION	Prevent large particulates and fluid exhaled by wearer from spreading
PERFORMANCE REQUIREMENTS	Bacterial filtration efficiency $\geq 98\%$, Particle Filtration Efficiency $\geq 98\%$ @ 0.1 micron, Differential pressure < 49 Pa, Fluid resistance ≥ 120 mmHg
FEATURES	Splash resistance, high breathability and comfort, hypoallergenic properties, tie-on design, no fit test required
SIZE	Standard - one size fits all
RECOMMENDED USE	Hospital - clinical use when risk of contact with airborne infections
REUSABLE	NO
FLUID RESISTANCE	Yes
AEROSOL PROTECTION	No
LATEX FREE	Yes
CERTIFICATIONS	CE, FDA