

Supporting health care around the globe.



Medical Surgical Masks (Type II)



Surgical masks are a barrier to splashes, droplets, and spit, and are not designed to seal tight against the face.

PRODUCT NAME	Surgical Masks (EN14683 - Type II)
PRODUCT TYPE	Fluid and droplet protection
PRODUCT DESCRIPTION	Prevent large particulates and minimal fluid exhaled by wearer from spreading
PERFORMANCE REQUIREMENTS	Bacterial filtration efficiency $\geq 98\%$, Particle Filtration Efficiency $\geq 98\%$ @ 0.1 micron, Differential pressure < 29 Pa, Fluid resistance ≥ 80 mmHg
FEATURES	High breathability and comfort, hypoallergenic properties, tie-on design, no fit test required
SIZE	Standard - one size fits all
RECOMMENDED USE	Hospitals - low risk clinical applications, not involving bloodborne pathogens or other bodily fluids
REUSABLE	No
FLUID RESISTANCE	No - Minimal
AEROSOL PROTECTION	No
LATEX FREE	Yes
CERTIFICATIONS	CE, FDA